Minimotion Set

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Match-U is a registered trademark of Match-U B.V. The Netherlands

WWW.MATCH-U.COM

NEDERLANDS

MINIMOTION

ENGLISH
DEUTSCH
FRANÇAIS
ESPAÑOL
NEDERLANDS

Match-U is a registered trademark of Match-U B.V. The Netherlands

WWW.MATCH-U.COM

Minimotion Set

Information for patients, physiotherapists and doctors.

Minimotion Set includes:

- Two Minimotion-holders with short band.
- One Minimotion-holder with medium band.
- One Minimotion-holder with large band.
- One light bar (28 cm).
- Four flexible bars with weights.

Attach the Minimotion around the body/device by using the straps. The confirmation should be firm but pleasant. Using the Minimotion can be practiced in three directions:

- forward / backward (holder on the forehead)
- sideways / (holder above the ears)
- rotation / (holder on the forehead)

Physiotherapists use the Minimotion for:

- Rehabilitation after surgery
- Rehabilitation after trauma
- instability
- Attitude weakness

Minimotion is not suitable for training:

- Within 48 hours after surgery or trauma.
- Operations where two vertebrae are fixed to each other (Fusion)
- Some tumors
- Fresh fractures

Side effects:

Side effects involving the body as a whole have not been reported to date. Condition is that the Minimotion handled properly and that the specified intensity is used.

The training possibilities:

Relaxation
During this kind of training there are very small forces generated.

Stabilization

During this form of training, small forces are generated. Similar to forces that occur during walking.

Muscle Strengthening

This form of training generates moderate to intense forces. Forces similar to isometric contractions and training with resistance bands.

Anticipation and speed training

During this kind of training you generate maximum forces. This form of training is only suitable for athletes and highly dutiable persons.

Choose the right intensity: www.match-u.com

If in doubt consult your doctor.

Instructions:

- For safety we recommend to keep a minimum distance of 1 meter from the Minimotion-user and other persons or objects. The flexible rod can cause damage when this minimum distance is not observed.
- Avoid sudden movements with the Minimotion. Unexpected movements of the Minimotion can injure/damage other people/objects.
- The flexibar in the holder should optimally be restrained. After tightening the screw the flexibar can not slide up or down anymore.
- The flexibar may not be on the other side of the holder stabbing.
- Children should only train with the Minimotion under adult supervision.
- Match-U B.V. is in no way responsible for any damage caused by the use of the Minimotion.

Washing instructions:

The various components of the Minimotion can be cleaned with a damp cloth, possibly in combination with a mild detergent.

	handwash only, maximum temperature 40°C handle with care
	do not use chlorine-based bleach
	do not tumble dry
	do not dry clean

Match-U is a registered trademark of Match-U B.V. The Netherlands

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