## Minimotion Set SMALL

NEDERLANDS

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WWW.MATCH-U.COM

Minimotion

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Information for patients, physiotherapists and doctors.

**Minimotion Set includes:**

- One Minimotion-holders with short band.
- One light bar (28 cm).
- Three flexible bars with weights.
- One clip.

Attach the Minimotion around the body/device by using the straps.

The confirmation should be firm but pleasant.

Using the Minimotion can be practiced in three directions:

- forward / backward (holder on the forehead)
- sideways / (holder above the ears)
- rotation / (holder on the forehead)

**Physiotherapists use the Minimotion for:**

- Rehabilitation after surgery
- Rehabilitation after trauma
- instability
- Attitude weakness

**Minimotion is not suitable for training:**

- Within 48 hours after surgery or trauma.
- Operations where two vertebrae are fixed to each other (Fusion)
- Some tumors
- Fresh fractures

**Side effects:**

Side effects involving the body as a whole have not been reported to date.

Condition is that the Minimotion handled properly and that the specified intensity is used.

**The training possibilities:**

Relaxation  
During this kind of training there are very small forces generated.

Stabilization  
During this form of training, small forces are generated.  
Similar to forces that occur during walking.

Muscle Strengthening  
This form of training generates moderate to intense forces.  
Forces similar to isometric contractions and training with resistance bands.

Anticipation and speed training  
During this kind of training you generate maximum forces.  
This form of training is only suitable for athletes and highly dutiable persons.

Choose the right intensity: [www.match-u.com](http://www.match-u.com)

*If in doubt consult your doctor.*

**Instructions:**

- For safety we recommend to keep a minimum distance of 1 meter from the Minimotion-user and other persons or objects. The flexible rod can cause damage when this minimum distance is not observed.
- Avoid sudden movements with the Minimotion. Unexpected movements of the Minimotion can injure/damage other people/objects.
- The flexibar in the holder should optimally restrained. After tightening the screw the flexibar can not slide up or down anymore.
- The flexibar may not be on the other side of the holder stabbing.
- Children should only train with the Minimotion under adult super-vision.
- Match-U B.V. is in no way responsible for any damage caused by the use of the Minimotion.

**Washing instructions:**

The various components of the Minimotion can be cleaned with a damp cloth, possibly in combination with a mild detergent.

	handwash only maximum temperature 40 C handle with care
	do not use chlorine-based bleach
	do not tumble dry
	do not dry clean

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